

INTEOS[®] - Humerus Plates 3.5



innovative distal Humerus Plating System



Advantages of the INTEOS® Humerus Plates 3.5 distal

The INTEOS® Humerus Plates 3.5 distal of Hofer-medical are characterized by its easy implantation technique as a result of pre-angled holes. **Fixed-angle & multidirectional** implantation techniques are as standard.

- + Angle stable and polyaxial
- + Anatomic precise shape
- + AO-approach (90°) as well as alternative-approach (180°) available
- + Chamfered auto-compression holes for generation of important compression
- + Screw heads are countersunk in the plate
- + Numerous possibilities for optimal temporary fragment fixation through Kirschner wire holes
- + Various plate lengths for maximum stable fracture supply
- + Milled groove on the bone side of the plate for medial and lateral optimal positioning of the plate along the supracondylar bone edge to find the correct position for the plate



INTEOS® Humerus distal dorsolateral 3.0 & 3.5 as av



Titanium	Number of Holes	Side	Length	PU
770-140-000-011	10 + 1 long hole	left	79 mm	1
770-140-000-014	13 + 1 long hole	left	103 mm	1
770-140-000-016	15 + 1 long hole	left	119 mm	1
770-140-000-018	17 + 1 long hole	left	135 mm	1
770-140-001-011	10 + 1 long hole	right	79 mm	1
770-140-001-014	13 + 1 long hole	right	103 mm	1
770-140-000-016	15 + 1 long hole	right	119 mm	1
770-140-000-018	17 + 1 long hole	right	135 mm	1

INTEOS® Humerus distal lateral 3.5 as av



Titanium	Number of Holes	Side	Length	PU
770-141-000-006	5 + 1 long hole	left	74 mm	1
770-141-000-008	7 + 1 long hole	left	90 mm	1
770-141-000-010	9 + 1 long hole	left	106 mm	1
770-141-000-012	11 + 1 long hole	left	122 mm	1
770-141-001-006	5 + 1 long hole	right	74 mm	1
770-141-001-008	7 + 1 long hole	right	90 mm	1
770-141-001-010	9 + 1 long hole	right	106 mm	1
770-141-001-012	11 + 1 long hole	right	122 mm	1

INTEOS® Humerus distal medial 3.5 as av



Titanium	Number of Holes	Side	Length	PU
770-142-000-008	7 + 1 long hole	left	88 mm	1
770-142-000-011	10+ 1 long hole	left	112 mm	1
770-142-000-012	10 + 2 long hole	left	129 mm	1
770-142-000-015	13 + 2 long hole	left	153 mm	1
770-142-001-008	7 + 1 long hole	right	88 mm	1
770-142-001-011	10+ 1 long hole	right	112 mm	1
770-142-001-012	10 + 2 long hole	right	129 mm	1
770-142-001-015	13 + 2 long hole	right	153 mm	1

INTEOS® Humerus distal medial 3.5 as av ext.



Titanium	Number of Holes	Side	Length	PU
770-143-000-009	8 + 1 long hole	left	94 mm	1
770-143-000-012	11 + 1 long hole	left	118 mm	1
770-143-000-013	11 + 2 long hole	left	135 mm	1
770-143-000-016	14 + 2 long hole	left	159 mm	1
770-143-001-009	8 + 1 long hole	right	94 mm	1
770-143-001-012	11 + 1 long hole	right	118 mm	1
770-143-001-013	11 + 2 long hole	right	135 mm	1
770-143-001-016	14 + 2 long hole	right	159 mm	1

Plates configuration and design in four models:
DL dorsolateral/ **L** lateral/ **M** medial / **M ext** medial extended



AUSTRIA

HOFER GMBH & CO KG
Jahnstrasse 10-12
A 8280 Fürstenfeld

TEL: +43 3382 53388
FAX: +43 3382 53093

office@hofer-medical.com
www.hofer-medical.com

FN21826y
UID ATU30764704



GERMANY

HOFER-MEDICAL GMBH
Ruhrtalstraße 52-60
D 45239 Essen

TEL: 0800-700 80 22 (innerhalb BRD)
FAX: +49 201 56599 1069

office@hofer-medical.de
www.hofer-medical.com

HRB24645
UID DE290596924



ITALY

HOFER MEDICAL ITALIA S.R.L.
Raingasse 26
39100 Bozen

TEL: +39 393 8950481
FAX: +39 06 81151148

office@hofer-medical.it
www.hofer-medical.it

RM-1454597
UID IT13534721009



POLAND

HOFER GMBH & CO KG SPÓŁKA
KOMANDYTOWA ODDZIAŁ W POLSCE
Radkowska 9
PL 57-402 Nowa Ruda

TEL: +48 603 224112
Fax: +48 748 712273

office@hofer-medical.pl
www.hofer-medical.com

REGON 141527293
UID PL1070011771

creating partnership

For us, a partnership is the result of a long and happy relationship in all areas of our work. Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your HOFER medical solutions Team

vCard

