


INTEOS[®] - Fibula Plate 3.0 & 3.5

innovative polyaxial distal Fibula Plating System



Advantages of the INTEOS® Fibula Plating System

The INTEOS® Fibula plate design showcases the latest insights into repositioning and stabilization of fractured fibulae, based on the **tried-and-tested fixed angle and polyaxial plate-screw-connection** by Hofer-medical.



Product features

- the plate design is based on over 100 precisely measured fibula bones
- precise anatomical plate prototype for best results
- multiple distal fixed-angle and polyaxial holes for optimized screw placement
- optional non-fixed-angle syndesmosis screw can be positioned in every hole on the plate
- minimized protrusion of screw head in order to avoid irritation and discomfort

Screw features

- compatible INTEOS® 3.0 & 3.5 mm locking screw
- Multi-directional, rounded screw head with locking function (fixed-angle & polyaxial)
- non-invasive screw tip
- „stick fit“ - the stable INTEOS® connection between screw and screwdriver blade
- ForceDrive No 8 & 10 - screwdriver technology

INTEOS® - the multidirectional fixed-angle screw system standard for small fragment screws



Stable anatomic design for safe preservation of repositioning



Accurate contour shape of the plate, based on over 100 measured fibulae



Innovative pre-bending technology of the plate in order to achieve best possible plate - bone alignment



AUSTRIA

HOFER GMBH & CO KG
Jahnstrasse 10-12
A 8280 Fürstenfeld

TEL: +43(0)3382/53388
FAX: +43(0)3382/53093

office@hofer-medical.com
www.hofer-medical.com

FN21826y
UID ATU30764704



GERMANY

HOFER-MEDICAL GMBH
Ruhrtalstraße 52-60
D 45239 Essen

TEL: 0800-700 80 22 (innerhalb BRD)
FAX: +49 201 56599 1069

office@hofer-medical.de
www.hofer-medical.com

HRB24645
UID DE290596924



ITALY

HOFER MEDICAL ITALIA S.R.L.
Via della Giuliana, 66
I 00195 Rome

Tel: +39 393 8950481
FAX: +39 06 81151148

office@hofer-medical.it
www.hofer-medical.it

RM-1454597
UID IT13534721009



POLAND

HOFER GMBH & CO KG SPÓŁKA
KOMANDYTOWA ODDZIAŁ W POLSCE
Radkowska 9
PL 57-402 Nowa Ruda

Tel.Kom.: +48603224112
Fax: +48 748712273

office@hofer-medical.pl
www.hofer-medical.com

REGON 141527293
UID PL1070011771

creating partnership

For us, a partnership is the result of a long and happy relationship in all areas of our work. Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your HOFER medical solutions Team

vCard

