


INTEOS® - Humerus Plate 3.5

innovative proximal Humerus Plating System



Advantages of the INTEOS® Humerus Plate 3.5 prox.

The INTEOS® Humerus Plate 3.5 prox. of Hofer-medical is characterized by its easy implantation technique as a result of pre-angled holes. **Fixed-angle & multidirectional** implantation techniques are as standard.

The INTEOS® Humerus Plate 3.5 prox. fulfills the high **ESS** requirements of Hofer-medical.

INTEOS® - Technology



Plate features

- Anatomically shaped plating system for the right and left humerus bone
- Orientation of the plate position nearby the sulcus bicipitalis for easy re-alignment of the fracture
- Pre-angulated holes with guided insertion device and optimized screw orientation
- Screw pairs crossed over one another ensure simultaneous higher resistance against misplacement of the head fragment and minimized risk of the screw moving into the joint
- Possibility of individual screw placement without the help of the guided insertion device but still with high angle flexibility
- Slimmer distal end of the plate for lower alteration of the deltoid muscle

Screw features

- Rounded screw head
- Non-invasive screw tip
- „Stick fit“ - the stable connection between screw driver and screw
- ForceDrive IO

Fixed-angle and multidirectional functionality as a result of the INTEOS® Small Fragment System standard



Exact anatomical fit with mandatory definition of the side



Increased stability of the system due to improved placement of the screws



Innovative screw placement for safe support of the joint



AUSTRIA

HOFER GMBH & CO KG
Jahnstrasse 10-12
A 8280 Fürstenfeld

TEL: +43 3382 53388
FAX: +43 3382 53093

office@hofer-medical.com
www.hofer-medical.com

FN21826y
UID ATU30764704



GERMANY

HOFER-MEDICAL GMBH
Ruhrtalstraße 52-60
D 45239 Essen

TEL: 0800-700 80 22 (innerhalb BRD)
FAX: +49 201 56599 1069

office@hofer-medical.de
www.hofer-medical.com

HRB24645
UID DE290596924



ITALY

HOFER MEDICAL ITALIA S.R.L.
Via della Giuliana, 66
I 00195 Rome

TEL: +39 393 8950481
FAX: +39 06 81151148

office@hofer-medical.it
www.hofer-medical.it

RM-1454597
UID IT13534721009



POLAND

HOFER GMBH & CO KG SPÓŁKA
KOMANDYTOWA ODDZIAŁ W POLSCE
Radkowska 9
PL 57-402 Nowa Ruda

TEL: +48 603 224112
Fax: +48 748 712273

office@hofer-medical.pl
www.hofer-medical.com

REGON 141527293
UID PL1070011771

creating partnership

For us, a partnership is the result of a long and happy relationship in all areas of our work. Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your HOFER medical solutions Team

vCard

